

PROVOX[®]

**Product Catalog
2017/2018**



Atos Atos Medical
Your voice

Atos Medical - Who we are

At Atos Medical, we are committed to giving a voice to people who breathe through a stoma, with design solutions and technologies built on decades of experience and a deep understanding of our users.

Since the introduction of the first Provox voice prosthesis in 1990, Atos Medical has become the world leader in Laryngectomy Care.

Our commitment and dedicated specialization in this field has allowed us to gain a deep understanding of the needs of our customers. This understanding drives our innovation process and we take pride in continuously bringing innovative solutions to our users.

Over the years we have served the needs of around 100,000 people with laryngectomy in over 70 countries. Our comprehensive Laryngectomy Care portfolio includes a complete range of Provox voice and pulmonary rehabilitation devices. This broad range of high quality Provox products allows for a personalized combination of solutions for each user.

Atos Medical understands that being the leader in this field comprises more than product development. Supporting clinical research and educating both professionals and users is an integral part of our business. Our products are featured in over 150 peer-reviewed publications. We support and conduct hundreds of conferences and educational programs every year.

We are proud of our history and our accomplishments and we continue to be excited by learning from our users; serving them with the best products, services, and education in Laryngectomy Care.

Patents and Trademarks:

Provox®, Provox® ActiValve®, Provox® FreeHands HME®, Provox® XtraBase® and Provox® TrachPhone® are registered trademarks owned by Atos Medical AB, Sweden.

Provox® Vega™, Provox® SmartInserter™, Provox® NID™, Provox® XtraMoist™ and XtraFlow™, Provox® Micron HME™, Provox® Luna™, Provox® HME Cap™, Provox® XtraFlange™, Provox® ShowerAid™, Provox® LaryTube™, Provox® LaryClip™ and Provox® LaryButton™, Provox® FlexiDerm™, Provox® OptiDerm™, Provox® FreeHands FlexiVoice™, Provox® FreeHands Support™, Provox® StabiliBase™ are trademarks of Atos Medical AB, Sweden.

For information about protective rights (e.g. patents), please refer to our webpage www.atosmedical.com/patents.

In this catalog you will find products that give the ability to speak, breathe well and take care of the skin.

Contents

4. Speaking

10. Breathing

17. Hands-free

21. Care

24. Laryngectomy illustrations

25. Selection of references

26. Contact

Speaking



Provox® Vega™

A voice prosthesis designed for good voice quality, low speaking effort, easy device maintenance and comfort.

- Less strained voice
- Easier to understand
- Produced with less effort
- Better voice quality



The Provox Vega voice prosthesis comes with an Insertion System that gives more control and greater flexibility.

- Adapts to patient anatomy and clinician's preferred insertion method
- Creates a perfect fold of the esophageal flange
- Intuitive and easy to use



Provox Vega 17Fr*	Size	Ref.no	Provox Vega 20Fr*	Size	Ref.no
Provox Vega 17Fr	4 mm	8270	Provox Vega 20Fr	4 mm	8276
Provox Vega 17Fr	6 mm	8271	Provox Vega 20Fr	6 mm	8277
Provox Vega 17Fr	8 mm	8272	Provox Vega 20Fr	8 mm	8278
Provox Vega 17Fr	10 mm	8273	Provox Vega 20Fr	10 mm	8279
Provox Vega 17Fr	12.5 mm	8274	Provox Vega 20Fr	12.5 mm	8280
Provox Vega 17Fr	15 mm	8275	Provox Vega 20Fr	15 mm	8281

Provox Vega 22.5Fr*	Size	Ref.no	Provox Capsule	Ref.no
Provox Vega 22.5Fr	4 mm	8282	Provox Capsule 17Fr	7795
Provox Vega 22.5Fr	6 mm	8283	Provox Capsule 20Fr	7796
Provox Vega 22.5Fr	8 mm	8284	Provox Capsule 22.5Fr	7797
Provox Vega 22.5Fr	10 mm	8285		
Provox Vega 22.5Fr	12.5 mm	8286		
Provox Vega 22.5Fr	15 mm	8287		

*Sterile, for single use.

* Includes: 1 Provox Vega voice prosthesis with Insertion System (sterile) and 1 Provox Brush (non-sterile).

Provox® Vega™ XtraSeal™

A solution to manage leakage around.

- The enlarged esophageal flange is thin and flexible which allows it to conform to the tissue around the TE puncture
- The user can continue to speak, eat and drink while the leakage around is managed



Provox Vega XtraSeal	Size	Ref.no	Provox Vega XtraSeal	Size	Ref.no	Provox Vega XtraSeal	Size	Ref.no
Provox Vega XtraSeal 17Fr	4 mm	8288	Provox Vega XtraSeal 20Fr	4 mm	8294	Provox Vega XtraSeal 22.5Fr	4 mm	8300
Provox Vega XtraSeal 17Fr	6 mm	8289	Provox Vega XtraSeal 20Fr	6 mm	8295	Provox Vega XtraSeal 22.5Fr	6 mm	8301
Provox Vega XtraSeal 17Fr	8 mm	8290	Provox Vega XtraSeal 20Fr	8 mm	8296	Provox Vega XtraSeal 22.5Fr	8 mm	8302
Provox Vega XtraSeal 17Fr	10 mm	8291	Provox Vega XtraSeal 20Fr	10 mm	8297	Provox Vega XtraSeal 22.5Fr	10 mm	8303
Provox Vega XtraSeal 17Fr	12.5 mm	8292	Provox Vega XtraSeal 20Fr	12.5 mm	8298	Provox Vega XtraSeal 22.5Fr	12.5 mm	8304
Provox Vega XtraSeal 17Fr	15 mm	8293	Provox Vega XtraSeal 20Fr	15 mm	8299	Provox Vega XtraSeal 22.5Fr	15 mm	8305

Sterile, for single use.

Includes: 1 Provox Vega XtraSeal prosthesis with Insertion System (sterile) and 1 Provox Brush (non-sterile).

Note: The shaft of the Provox XtraSeal is ca 1 mm shorter than the size indicated due to the enlarged esophageal flange.

Provox® ActiValve®

The voice prosthesis with extended device lifetime.

- Increased device life – reduces the need of frequent replacements in most patients who experience short device life
- Available in three strengths – Light, Strong and XtraStrong



Provox ActiValve	Size	Ref.no	Provox ActiValve	Size	Ref.no	Provox ActiValve	Size	Ref.no
Provox ActiValve Light	4.5 mm	7150	Provox ActiValve Strong	4.5 mm	7160	Provox ActiValve XtraStrong	4.5 mm	7165
Provox ActiValve Light	6 mm	7151	Provox ActiValve Strong	6 mm	7161	Provox ActiValve XtraStrong	6 mm	7166
Provox ActiValve Light	8 mm	7152	Provox ActiValve Strong	8 mm	7162	Provox ActiValve XtraStrong	8 mm	7167
Provox ActiValve Light	10 mm	7153	Provox ActiValve Strong	10 mm	7163	Provox ActiValve XtraStrong	10 mm	7168
Provox ActiValve Light	12.5 mm	7154	Provox ActiValve Strong	12.5 mm	7164	Provox ActiValve XtraStrong	12.5 mm	7169
						Provox ActiValve Lubricant (replacement)		7149

Non sterile, for single use. Includes: 1 Provox ActiValve voice prosthesis, 1 single-use insertion tool, 2 Provox Brush, 1 Provox ActiValve Lubricant, 1 Provox Plug, 2 User Cards, 1 Emergency Card.

Provox® Vega™ Puncture Set

A disposable surgical set for primary and secondary placement of Provox Vega.

- All-in-one surgical set
- Creates a round puncture, leaving the surrounding tissue intact
- Comes pre-loaded with a Provox Vega voice prosthesis
- Disposable and sterile
- Enables a rapid procedure, minimum preparation and assembly are needed



Provox Vega Puncture Set	Size	Ref. no
Provox Vega Puncture Set 17 Fr	8 mm	8140
Provox Vega Puncture Set 17 Fr	10 mm	8141
Provox Vega Puncture Set 17 Fr	12.5 mm	8142
Provox Vega Puncture Set 17 Fr	15 mm	8143
Provox Vega Puncture Set 20 Fr	8 mm	8144

Provox Vega Puncture Set	Size	Ref. no
Provox Vega Puncture Set 20 Fr	10 mm	8145
Provox Vega Puncture Set 20 Fr	12.5 mm	8146
Provox Vega Puncture Set 22.5 Fr	8mm	8147
Provox Vega Puncture Set 22.5 Fr	10 mm	8148
Provox Vega Puncture Set 22.5 Fr	12.5 mm	8149

Includes: 1 Pharynx Protector (made of transparent thermoplastic), 1 Puncture Needle (made of surgical stainless steel), 1 Guidewire (made of pre-colored fluoroplastic), 1 Puncture Dilator with 1 preloaded Provox Vega voice prosthesis and 1 Provox Brush.

Provox® Brush

For cleaning of the voice prosthesis.

- Easy and efficient cleaning
- Can be bent for individual adaptation
- Can be used to insert Provox Plug and Provox Vega Plug

Provox Brush	Quantity	Size	Ref.no
Provox Brush	6 pcs	4.5-10 mm	7204
Provox Brush XL	6 pcs	12-18 mm	7225

Non sterile, reusable. For single patient use.
For voice prostheses in sizes 12.5 and 15 mm use Provox Brush XL.



Provox® Flush

For flushing after brushing.

- Can be used with air or drinking water
- May be bent for individual fitting

Provox Flush	Quantity	Ref.no
Provox Flush	1 pc	8109

Non sterile, reusable. For single patient use.



Provox® Plug & Provox® Vega™ Plug

A first-aid to temporarily stop leakage through the voice prosthesis.

Provox Plug & Provox Vega Plug	Quantity	Ref.no
Provox Plug	1 pc	7205
Provox Vega Plug, 17 Fr	1 pc	8119
Provox Vega Plug, 20 Fr	1 pc	8129
Provox Vega Plug, 22.5 Fr	1 pc	8139

Non sterile, reusable. For single patient use.
Includes Provox Brush.

Note: The Provox Vega Plug fits into Provox Vega prostheses and the Provox Plug is for Provox2, Provox (1) and Provox ActiValve.




Provox® XtraFlange™

A rapid solution to manage leakage around.

Provox XtraFlange	Quantity	Ref.no
Provox XtraFlange 22.5 Fr	1 pc	7275
Provox XtraFlange 20 Fr	1 pc	7276
Provox XtraFlange 17 Fr	1 pc	7277

Sterile, for single use.



Provox® GuideWire

Provox GuideWire enables retrograde insertion of all Provox indwelling prostheses. The tip contains a retainer slit that holds the safety strap of the prosthesis.

Provox GuideWire	Quantity	Ref.no
Provox GuideWire	1 pc	7215

Sterile, for single use.



Provox® Measure

Provox Measure is used to determine the length (corresponding to prosthesis size) of TE punctures. The correct prosthesis size (1-2 mm longer than the puncture) is important to minimize the risk of complications.

Provox Measure	Quantity	Ref.no
Provox Measure	1 pc	7270
Provox Measure Flanges	5 pcs	7271

Provox Measure: non sterile, reusable, sterilizable.
Provox Measure Flanges: non sterile, for single use.



Provox® Dilator

Provox Dilator is used to upsize the diameter of the TE puncture or to temporary block a TE puncture, irrespective of type and diameter of voice prosthesis used. The size can be increased from 15 Fr to 24 Fr in one device (steps: 18, 20, 22, 24 Fr).

Provox Dilator	Quantity	Ref.no
Provox Dilator	1 pc	7211
Provox Dilator 17	1 pc	7122
Provox Dilator 20	1 pc	7123

Non sterile, reusable, sterilizable. Two smaller dilators are also available; for dilating 17 Fr or 20 Fr punctures.



Breathing



Provox® XtraFlow™

This HME helps to make it easier to breathe. Use it during the day when being physically active.

- Low breathing resistance
- Reduces coughing and mucus production
- Easy occlusion of the stoma for voicing
- Low-profile design



Provox XtraHME	Quantity	Ref.no
Provox XtraFlow HME	30 pcs	7291

Provox® XtraMoist™

This HME makes the air more humid. Use it during the day, when taking it easy.

- Comfortable breathing
- Reduces coughing and mucus production
- Easy occlusion of the stoma for voicing
- Low-profile design



Provox XtraHME	Quantity	Ref.no
Provox XtraMoist HME	30 pcs	7290

Provox® Luna™

A night-time solution to help sleep comfortably, soothe the skin and improve lung health.

- Soft and smooth for night-time comfort
- Adhesive made of hydrogel to soothe and cool the skin
- Superior humidification compared to other HMEs*
- Low breathing resistance for easy night-time breathing
- Side openings to prevent occlusion while sleeping



Provox Luna	Quantity	Ref.no
Provox Luna HME	30 pcs	8013
Provox Luna Adhesive	15 pcs	8014

* Data on file

Provox® Micron HME™

Not only an HME, but also an effective bacteria, virus, and dust filter.

- A highly effective bacterial and viral filter (≥99%)
- Reduces coughing and mucus production

Please note that the Provox Micron filtration capacity may vary depending on the size, shape and speed of airborne particles.



Provox Micron HME	Quantity	Ref.no
Provox MicronHME	5 pcs	7247
Provox MicronHME	30 pcs	7248

Non sterile, for single use.

Provox® StabiliBase™

Gives superior seal and support for the stoma when speaking hands-free.

- For speaking hands-free
- For speaking much with finger occlusion
- Sticks well to the skin
- For a deep stoma

Provox StabiliBase	Quantity	Ref.no
Provox StabiliBase	15 pcs	7289
Provox StabiliBase OptiDerm	15 pcs	7318

Non sterile, for single use.



Provox® XtraBase®

For speaking with finger occlusion or hands-free.

- For speaking with finger occlusion or hands-free
- For a flat or somewhat deeper stoma
- Sticks well to the skin

Provox XtraBase	Quantity	Ref.no
Provox XtraBase	20 pcs	7265

Non sterile, for single use.




Provox® FlexiDerm™

Gives a good seal around the stoma during the day.

- For speaking with finger occlusion
- Sticks well to the skin
- For a flat stoma

Provox FlexiDerm	Quantity	Ref.no
Provox FlexiDerm Round	20 pcs	7253
Provox FlexiDerm Oval	20 pcs	7254
Provox FlexiDerm Plus	20 pcs	7331

Non sterile, for single use.



Provox® OptiDerm™

An adhesive for the day if the skin feels sore or irritated.

- For sore or irritated skin
- Use after surgery
- Use after radiotherapy

Provox OptiDerm	Quantity	Ref.no
Provox OptiDerm Round	20 pcs	7255
Provox OptiDerm Oval	20 pcs	7256
Provox OptiDerm Plus	20 pcs	7332

Non sterile, for single use.



Provox® LaryButton™

For maintaining the opening of the stoma.

Provox LaryButton	OD	L	Ref.no
Provox LaryButton 12/8	12	8	7671
Provox LaryButton 14/8	14	8	7672
Provox LaryButton 16/8	16	8	7673
Provox LaryButton 18/8	18	8	7674

Provox LaryButton	OD	L	Ref.no
Provox LaryButton 12/18	12	18	7685
Provox LaryButton 14/18	14	18	7686
Provox LaryButton 16/18	16	18	7687
Provox LaryButton 18/18	18	18	7688

Non sterile, for single patient use.



Provox® LaryTube™

For maintaining the opening of the tracheostoma.
Available in Standard, Fenestrated and with Ring
(to be secured with Provox adhesive).

ID = inner diameter
OD = outer diameter
L = length (in mm)

Provox® LaryTube™ Standard

Standard	ID	OD	L	Ref.no	Standard	ID	OD	L	Ref.no
Provox LaryTube 8/27	9.5	12.0	27	7601	Provox LaryTube 10/27	12.0	15.0	27	7609
Provox LaryTube 8/36	9.5	12.0	36	7602	Provox LaryTube 10/36	12.0	15.0	36	7610
Provox LaryTube 8/55	9.5	12.0	55	7603	Provox LaryTube 10/55	12.0	15.0	55	7611
Provox LaryTube 9/27	10.5	13.5	27	7605	Provox LaryTube 12/27	13.5	17.0	27	7613
Provox LaryTube 9/36	10.5	13.5	36	7606	Provox LaryTube 12/36	13.5	17.0	36	7614
Provox LaryTube 9/55	10.5	13.5	55	7607	Provox LaryTube 12/55	13.5	17.0	55	7615

Includes: 5 XtraFlow Cassettes (7298).
Non sterile.



Provox® LaryTube™ Fenestrated

Fenestrated	ID	OD	L	Ref.no	Fenestrated	ID	OD	L	Ref.no
Provox LaryTube 8/36	9.5	12.0	36	7637	Provox LaryTube 10/36	12.0	15.0	36	7643
Provox LaryTube 8/55	9.5	12.0	55	7638	Provox LaryTube 10/55	12.0	15.0	55	7644
Provox LaryTube 9/36	10.5	13.5	36	7640	Provox LaryTube 12/36	13.5	17.0	36	7646
Provox LaryTube 9/55	10.5	13.5	55	7641	Provox LaryTube 12/55	13.5	17.0	55	7647

Includes: 5 XtraFlow Cassettes (7298) and
1 Provox Brush (7204).
Non sterile.



Provox® LaryTube™ with Ring

With Ring	ID	OD	L	Ref.no	With Ring	ID	OD	L	Ref.no
Provox LaryTube 8/36	9.5	12.0	36	7624	Provox LaryTube 10/36	12.0	15.0	36	7628
Provox LaryTube 8/55	9.5	12.0	55	7625	Provox LaryTube 10/55	12.0	15.0	55	7629
Provox LaryTube 9/36	10.5	13.5	36	7626	Provox LaryTube 12/36	13.5	17.0	36	7630
Provox LaryTube 9/55	10.5	13.5	55	7627	Provox LaryTube 12/55	13.5	17.0	55	7631

Includes: 5 XtraFlow Cassettes (7298),
1 FlexiDerm Oval (7254).
Non sterile.



Provox® TubeHolder

The Provox TubeHolder has been developed exclusively for use with Provox LaryTube and Provox LaryButton. It is suitable for users with sensitive skin around the tracheostoma. The integrated clip connectors allow for optimal fit to the wings of the Provox LaryTube and LaryButton.

Provox TubeHolder	Quantity	Ref.no
Provox TubeHolder	1 pc	7668

Non sterile. For single use.



Provox® LaryClip™

The Provox LaryClip is an alternative to Provox TubeHolder.

Provox LaryClip	Ref.no
1 set: 8 pcs LaryClip, 40 pcs LaryClip Base	7669

Non sterile, for single use.



Provox® HME Cap

Provox HME Cap is used with the FreeHands HMEs. The Cap can be occluded manually for speech.

Provox HME Cap	Quantity	Ref.no
Provox HME Cap matt titanium	1 pc	7730

Reusable, non sterile.



Provox® TubeBrush

The Provox TubeBrush is used to clean the Provox LaryTube and Provox LaryButton.

Provox TubeBrush	Quantity	Size	Ref.no
Provox TubeBrush	6 pcs	8 mm	7660
Provox TubeBrush	6 pcs	12 mm	7661

Non sterile, reusable. For single patient use.



Provox® BasePlate Adaptor

Provox BasePlate Adaptor makes it possible to attach a device, e.g. an HME, with an ISO 15 mm standard connector to a Provox attachment.

Provox BasePlate Adaptor	Quantity	Ref.no
Provox BasePlate Adaptor	1 pc	7263

Non sterile, for single patient use.



Provox® HME Cassette Adaptor

Provox HME Cassette Adaptor enables use of a Provox HME Cassette to any tracheostomy tube with a 15 mm ISO connector.

Provox HME Cassette Adaptor	Quantity	Ref.no
Provox HME Cassette Adaptor	1 pc	7246

Non sterile, for single patient use.



Provox® FenestrationPunch

The Provox FenestrationPunch is used to make small fenestrations in the Provox LaryTube.

Provox FenestrationPunch	Quantity	Ref.no
Provox FenestrationPunch	1 pc	7654

Non sterile, reusable.



Provox® LaryButton™ Sizer Kit and Provox® LaryTube™ Sizer Kit

To allow for proper selection and sizing, the Provox LaryButton Sizer Kit and Provox LaryTube Sizer Kit contains sterilizable and reusable tubes for individual fitting. The products can be secured with Provox LaryClip and Provox TubeHolder.

Provox LaryButton Sizer Kit	Ref.no
Provox LaryButton Sizer Kit	7690

Non sterile, reusable, sterilizable (Sizer Buttons and inner box). Outer tray can be cleaned in instrument washer.

Provox LaryTube Sizer Kit	Ref.no
Provox LaryTube Sizer Kit	7648

Non sterile, reusable, sterilizable (Sizer Tubes and inner box). Outer tray can be cleaned in instrument washer.



Hands-free



Provox® FreeHands FlexiVoice™

Gives the freedom to speak hands-free.
Available in four membranes;
Light, Medium, Strong, XtraStrong.

- Speak hands-free speech or with finger occlusion when wanted
- Use different membranes for different situations
- Lock for physical activity
- Fast set up
- Reduces mucus production and coughing



Provox FreeHands FlexiVoice	Ref.no	Provox FreeHands FlexiVoice	Ref.no
Provox FreeHands FlexiVoice Set Plus: 1 light, 1 medium, 1 strong, 2 Arches, 1 Removal Aid, 30 HME Cassettes	7757	Provox FreeHands FlexiVoice Light: 1 light, 1 Arch, 1 Removal Aid	8161
Provox FreeHands FlexiVoice Set: 1 light, 1 medium, 1 strong, 2 Arches, 1 Removal Aid	7760	Provox FreeHands FlexiVoice Medium: 1 medium, 1 Arch, 1 Removal Aid	8162
		Provox FreeHands FlexiVoice Strong: 1 strong, 1 Arch, 1 Removal Aid	8163
		Provox FreeHands FlexiVoice XtraStrong: 1 XtraStrong, 1 Arch, 1 Removal Aid	8166

Provox® FreeHands Support™

To reduce stoma movement and improve voice quality when speaking hands-free with Provox FreeHands FlexiVoice.

Provox FreeHands Support	Quantity	Ref.no
Provox FreeHands Support Starter Set: 1 Support Flat, 1 Support Medium, 1 Support Deep	1 pc	8020
Provox FreeHands Support Flat	1 pc	8021
Provox FreeHands Support Medium	1 pc	8022
Provox FreeHands Support Deep	1 pc	8023
Provox FreeHands Support Adhesive	15 pcs	8024



Provox® FreeHands HME® Flow

For humidification when being active.

- To be combined with Provox FreeHands FlexiVoice when speaking hands-free

Provox FreeHands HME Flow	Quantity	Ref.no
Provox FreeHands HME Flow	30 pcs	8221

Non sterile, for single use.



Provox® FreeHands HME® Moist

For extra humidification when relaxing.

- To be combined with Provox FreeHands FlexiVoice when speaking hands-free

Provox FreeHands HME Moist	Quantity	Ref.no
Provox FreeHands HME Moist	30 pcs	8220

Non sterile, for single use.



Provox® FreeHands FlexiVoice™ Arch

To prevent fabric from occluding FlexiVoice.

Provox FreeHands FlexiVoice Arch	Quantity	Ref.no
Provox FreeHands FlexiVoice Arch	5 pcs	8165



Provox® StabiliBase™

Gives superior seal and support for the stoma when speaking hands-free.

- For speaking hands-free
- For speaking much with finger occlusion
- Sticks well to the skin
- For a deep stoma

Provox StabiliBase	Quantity	Ref.no
Provox StabiliBase	15 pcs	7289
Provox StabiliBase OptiDerm	15 pcs	7318

Non sterile, for single use.



Provox® XtraBase®

For speaking with finger occlusion or hands-free.

- For speaking with finger occlusion or hands-free
- For a flat or somewhat deeper stoma
- Sticks well to the skin

Provox XtraBase	Quantity	Ref.no
Provox XtraBase	20 pcs	7265

Non sterile, for single use.



Care



Provox® Adhesive Strip™

Provides additional adhesive seal. Protects Provox Luna when showering.

Provox Adhesive Strip	Quantity	Ref.no
Provox Adhesive Strip	10 pcs	8015



Provox® Luna™ ShowerAid

Protects the stoma from water when showering with Provox Luna.

Provox Luna ShowerAid	Quantity	Ref.no
Provox Luna ShowerAid	1 pc	8016



Provox® Adhesive Remover

To help removing the adhesive.

Provox Adhesive Remover	Quantity	Ref.no
Provox Adhesive Remover	50 pcs	8012

Non sterile, for single use.



Provox® Skin Barrier

To leave a protective layer on the skin.

Provox Skin Barrier	Quantity	Ref.no
Provox Skin Barrier	50 pcs	8011

Non sterile, for single use.



Provox® Cleaning Towel

To clean the skin before putting on an adhesive.

Provox Cleaning Towel	Quantity	Ref.no
Provox Cleaning Towel	10 pcs	7244

(20 resealable pouches of 10).
Non sterile, for single use.



Skin Tac™

Improves adhesion of appliances and tape to the skin.

SkinTac	Quantity	Ref.no
Skin Tac Adhesive Barrier Wipes	50 pcs	MS407W



Removal Aid

To remove the HME from the FlexiVoice speaking valve.

Removal Aid	Quantity	Ref.no
Removal Aid	2 pcs	7745



Provox® Silicone Glue

Improves the adhesion between the skin and the adhesive.

Provox Silicone Glue	Quantity	Ref.no
Provox Silicone Glue	1 pc	7720

Bottle with brush cap, non sterile.
For multiple application.



Provox® ShowerAid


Protects the stoma from water when showering.

Provox ShowerAid	Quantity	Ref.no
Provox ShowerAid	1 pc	7260


Non sterile, for single patient use.
Includes: 1 Provox FlexiDerm Oval Adhesive (7254).




Laryngectomy illustrations




Before total laryngectomy



After total laryngectomy



Speaking with the Provox voice prosthesis:
Manual occlusion of the Provox XtraHME
attached with Provox Adhesive.



Speaking with the Provox voice prosthesis:
Handsfree with Provox FreeHands FlexiVoice
attached with Provox Adhesive.

Selection of references

Provox® Voice Prostheses

- Timmermans AJ, Harmsen HJ, Bus-Spoor C, Buijssen KJ, van As-Brooks C, de Goffau MC, Tonk RH, van den Brekel MW, Hilgers FJ, van der Laan BF. Biofilm formation on the Provox ActiValve: Composition and ingrowth analyzed by Illumina paired-end RNA sequencing, fluorescence in situ hybridization, and confocal laser scanning microscopy. *Head Neck*. 2016 Apr;38 Suppl 1:E432-40.
- Kress P, Schäfer P, Schwerdtfeger FP, Rösler S. Are modern voice prostheses better? A lifetime comparison of 749 voice prostheses. *Eur Arch Otorhinolaryngol*. 2014 Jan;271(1):133-40.
- Hancock KL, Lawson NR, Ward EC. Device life of the Provox Vega voice prosthesis. *Eur Arch Otorhinolaryngol*. 2013 Mar;270(4):1447-53.
- Hilgers FJ, Lorenz KJ, Maier H et al. Development and (pre-) clinical assessment of a novel surgical tool for primary and secondary tracheoesophageal puncture with immediate voice prosthesis insertion, the Provox Vega Puncture Set. *Eur Arch Otorhinolaryngol*. 2013 Jan;270(1):255-62.
- Hancock K, Ward E, Lawson N, van As-Brooks CJ. A prospective, randomized comparative study of patient perceptions and preferences of two types of indwelling voice prostheses. *Int J Lang Commun Disord*. 2012 May-Jun;47(3):300-9.
- Ward EC, Hancock K, Lawson N, van As-Brooks CJ. Perceptual characteristics of speech production using the new indwelling Provox Vega voice prosthesis: A randomized controlled cross-over trial. *Head & Neck*. 2011 Jan;33(1):13-9.
- Hilgers FJ, Ackerstaff AH, Jacobi I, Balm AJ, Tan B, Van den Brekel MW. Prospective clinical phase II study of two new indwelling voice prostheses (Provox Vega 22.5 and 20 Fr) and a novel anterograde insertion device (Provox SmartInserter). *Laryngoscope*. 2010 Jun;120(6):1135-43.
- Hilgers FJ, Ackerstaff AH, van RM et al. Clinical phase I/feasibility study of the next generation indwelling Provox voice prosthesis (Provox Vega). *Acta Otolaryngologica*. 2010 April; 130(4):511-9.
- Soolsma J, van den Brekel MW, Ackerstaff AH, Balm AJ, Tan B, Hilgers FJ. Long-Term Results of Provox ActiValve, Solving the Problem of Frequent Candida- and "Underpressure"-Related Voice Prosthesis Replacements. *The Laryngoscope*. 2007 The American Laryngological, Rhinological and Otolological Society, Inc.
- Hancock K, Houghton B, Van As-Brooks CJ, Coman W. First clinical experience with a new nonindwelling voice prosthesis (Provox NID) for voice rehabilitation after total laryngectomy. *Acta Otolaryngologica*. 2005 Sep;125(9):981-90.
- Hilgers FJ, Ackerstaff AH, Balm AJ, Van den Brekel MW, Bing Tan I, Persson JO. A new problemsolving indwelling voice prosthesis, eliminating the need for frequent Candida- and "underpressure"-related replacements: Provox ActiValve. *Acta Otolaryngologica*. 2003 Oct;123(8):972-9.
- Free RH, Van der Mei HC, Elving GJ, Van Weissenbruch R, Albers FW, Busscher HJ. Influence of the Provox Flush, blowing and imitated coughing on voice prosthetic biofilms in vitro. *Acta Otolaryngologica*. 2003 May;123(4):547-51.
- Op de Coul BMR, Hilgers FJM, Balm AJM, Tan IB, van den Hoogen FJA, van Tinteren H. A decade of post-laryngectomy vocal rehabilitation in 318 patients: A single institution's experience with consistent application of Provox® indwelling voice prostheses. *Archives of Otolaryngology- Head & Neck Surgery*. 2000 Nov;126(11):1320-8.
- Balle VH, Rindø L, Thomsen JC. Primary speech restoration at laryngectomy by insertion of voice prosthesis - 10 years experience. *Acta Otolaryngologica* 2000; Suppl 543:244-245.
- Ackerstaff AH, Hilgers FJ, Meeuwis CA, van der Velden LA, van den Hoogen FJ, Marres HA, Vreeburg GC, Manni JJ. Multi-institutional assessment of the Provox®2 voice prosthesis. *Archives of Otolaryngology - Head & Neck Surgery*. 1999 Feb;125(2):167-73.
- Hilgers FJM, Ackerstaff AH, Balm AJM, Tan IB, Aaronson NK, Persson JO. Development and clinical evaluation of a second-generation voice prosthesis (Provox®2) designed for anterograde and retrograde insertion. *Acta Otolaryngologica*. 1997 Nov;117(6):889-96.
- Hilgers FJM, Schouwenburg PF. A new low-resistance, self-retaining prosthesis (Provox®) for voice rehabilitation after total laryngectomy. *Laryngoscope*. 1990 Nov;100(11):1202-7.

Provox® HME and Provox FreeHands HME

- Lansaat L, de Kleijn BJ, Hilgers FJ, van der Laan BF, van den Brekel MW. A prospective multicenter clinical feasibility study of a new automatic speaking valve for postlaryngectomy voice rehabilitation. *Eur Arch Otorhinolaryngol*. 2017 Feb;274(2):1005-1013.
- Macri GF, Bogaardt H, Parrilla C, Minni A, D'Alatri L, de Vincentiis M, Paludetti G. Patients' experiences with HMEs and attachments after total laryngectomy. *Clin Otolaryngol*. 2016 Dec;41(6):652-659.
- van den Boer C, Muller SH, van der Noort V, Olmos RA, Minni A, Parrilla C, Hilgers FJ, van den Brekel MW, van der Baan S. Effects of heat and moisture exchangers on tracheal mucociliary clearance in laryngectomized patients: a multi-center case-control study. *Eur Arch Otorhinolaryngol*. 2015 Nov;272(11):3439-50.
- Parrilla C, Minni A, Bogaardt H, Macri GF, Battista M, Roukos R, Pandolfini M, Ruoppolo G, Paludetti G, D'Alatri L, de Vincentiis M. Pulmonary Rehabilitation After Total Laryngectomy: A Multicenter Time-Series Clinical Trial Evaluating the Provox XtraHME in HME-Naïve Patients. *Ann Otol Rhinol Laryngol*. 2015 Sep;124(9):706-13.
- van den Boer C, Muller SH, Vincent AD, van den Brekel MW, Hilgers FJ. Ex vivo assessment and validation of water exchange performance of 23 Heat and Moisture Exchangers for laryngectomized patients. *Respir Care*. 2014 Aug;59(8):1161-71.
- van den Boer C, Muller SH, Vincent AD, Züchner K, van den Brekel MW, Hilgers FJ. Ex vivo water exchange performance and short-term clinical feasibility assessment of newly developed heat and moisture exchangers for pulmonary rehabilitation after total laryngectomy. *Eur Arch Otorhinolaryngol*. 2014 Feb;271(2):359-66.
- Herranz J, Espiño MA, Morado CO. Pulmonary rehabilitation after total laryngectomy: a randomized cross-over clinical trial comparing two different heat and moisture exchangers (HMEs). *Eur Arch Otorhinolaryngol*. 2013 Sep;270(9):2479-84.
- Brook I, Bogaardt H, van As-Brooks C. Long-term use of heat and moisture exchangers among laryngectomies: medical, social, and psychological patterns. *Ann Otol Rhinol Laryngol*. 2013 Jun;122(6):358-63.
- Hilgers FJ, Dirven R, Wouters Y, Jacobi I, Marres HA, van den Brekel MW. A multicenter, prospective, clinical trial evaluating a novel adhesive baseplate (Provox StabiliBase) for peristomal attachment of postlaryngectomy pulmonary and voice rehabilitation devices. *Laryngoscope*. 2012 Nov;122(11):2447-53.
- Scheenstra RJ, Muller SH, Hilgers FJ. Endotracheal temperature and humidity in laryngectomized patients in a warm and dry environment and the effect of a heat and moisture exchanger. *Head Neck*. 2011 Sep;33(9):1285-93.
- Scheenstra RJ, Muller SH, Vincent A, Ackerstaff AH, Jacobi I, Hilgers FJ. A new heat and moisture exchanger for laryngectomized patients: endotracheal temperature and humidity. *Respir Care*. 2011 May;56(5):604-11.
- Scheenstra RJ, Muller SH, Vincent A, Hilgers FJ. Heat and moisture exchange capacity of the upper respiratory tract and the effect of tracheotomy breathing on endotracheal climate. *Head Neck*. 2011 Jan;33(1):117-24.
- Scheenstra RJ, Muller SH, Vincent A, Sinaasappel M, Hilgers FJ. Influence of breathing resistance of heat and moisture exchangers on tracheal climate and breathing pattern in laryngectomized individuals. *Head & Neck*. 2010 Aug;32(8):1069-78.
- Scheenstra RJ, Muller SH, Vincent A, Ackerstaff AH, Jacobi I, Hilgers FJM. Short term endotracheal climate changes and clinical effects of a heat and moisture exchanger with an integrated electrostatic virus and bacterial filter developed for laryngectomized individuals. *Acta Otolaryngologica*. 2010 Jun;130(6):739-46.
- Bien S, Okla S, van As-Brooks CJ, Ackerstaff AH. The effect of a Heat and Moisture Exchanger (Provox HME) on pulmonary protection after total laryngectomy: a randomized controlled study. *European Archives of Otorhinolaryngology*. 2010 March;267(3):429-35.
- Lorenz KJ, Huverstuhl K, Maier H. [Finger-free speech with the Provox FreeHands HME Automatic Tracheostoma Valve system. Clinical long-term experience]. *HNO*. 2009 Nov;57(11):1090-8.

- Zuur JK, Muller SH, Vincent A, Sinaasappel M, de Jongh FH, Hilgers FJ. The influence of a heat and moisture exchanger on tracheal climate in a cold environment. *Medical & Biological Engineering & Computing*. 2009 Sep;31(7):852-7.
- Lorenz KJ, Maier H. [Pulmonary rehabilitation after total laryngectomy using a heat and moisture exchanger (HME)]. *Laryngorhinootologie*. 2009 Aug;88(8):513-22.
- Scheenstra RJ, Muller SH, Vincent A, Sinaasappel M, Zuur JK, Hilgers FJ. Endotracheal temperature and humidity measurements in laryngectomized patients: intra- and inter-patient variability. *Medical & Biological Engineering & Computing*. 2009 Jul;47(7):773-82.
- Zuur JK, Muller SH, Vincent A, Sinaasappel M, de Jongh FH, Hilgers FJ. Assessment of tracheal temperature and humidity in laryngectomized individuals and the influence of a heat and moisture exchanger on tracheal climate. *Head & Neck*. 2008 Aug;30(8):1072-1082.
- Hilgers FJ, Soolsma J, Ackerstaff AH, Balm AJ, Tan IB, van den Brekel MW. A thin tracheal silicone washer to solve periprosthetic leakage in laryngectomies: direct results and long-term clinical effects. *Laryngoscope*. 2008 Apr;118(4):640-5.
- Zuur JK, Muller SH, Sinaasappel M, Hart GA, Van Zandwijk N, Hilgers FJ. Influence of heat and moisture exchanger respiratory load on transcutaneous oxygenation in laryngectomized individuals: a randomized cross-over study. *Head & Neck*. 2007 Dec;29(12):1102-1110.
- Lorenz KJ, Ackerstaff AH, Hilgers FJM, Maier H. Hands-free speech after surgical voice rehabilitation with a Provox® voice prosthesis: experience with the Provox® Freehands™ HME tracheostoma valve system. *European Archives of Otorhinolaryngology*. 2007 Feb; 264(2):151-157.
- Hilgers FJ, Ackerstaff AH, van RM et al. Development and evaluation of a novel tracheostoma button and fixation system (Provox LaryButton and LaryClip adhesive) to facilitate hands-free tracheoesophageal speech. *Acta Otolaryngologica*. 2006 Dec;126(11):1218-1224.
- Hamade R, Hewlett N, & Scanlon, E. A quantitative and qualitative evaluation of an automatic occlusion device for tracheoesophageal speech: the Provox FreeHands HME. *Clin Linguist Phon*. 2006 Apr-May;20(2-3):187-93.
- Zuur JK, Muller SH, de Jongh FH, van Zandwijk N, Hilgers FJ. The physiological rationale of heat and moisture exchangers in post-laryngectomy pulmonary rehabilitation: a review. *European Archives of Otorhinolaryngology*. 2006 Jan;263(1):1-8.
- Tervonen H, Back L, Juvas A, Rasanen P, Makitie AA, Sintonen H, Roine RP, Vilkanen E, Aaltonen LM. Automatic speaking valve in speech rehabilitation for laryngectomized patients. *European Archives of Otorhinolaryngology*. 2005 Oct;262(10):816-20.
- Op de Coul BM, Ackerstaff AH, van As-Brooks CJ, van den Hoogen FJ, Meeuwis CA, Manni JJ, Hilgers FJ. Compliance, quality of life and quantitative voice quality aspects of hands-free speech. *Acta Otolaryngologica*. 2005 Jun;125(6):629-37.
- Ackerstaff AH, Fuller D, Irvin M, Maccracken E, Gaziano J, Stachowiak L. Multicenter study assessing effects of heat and moisture exchanger use on respiratory symptoms and voice quality in laryngectomized individuals. *Otolaryngology Head & Neck Surgery*. 2003 Dec;129(6):705-12.
- Hilgers FJ, Ackerstaff AH, Van As CJ, Balm AJ, Van den Brekel MW, Tan IB. Development and clinical assessment of a heat and moisture exchanger with a multi-magnet automatic tracheostoma valve (Provox FreeHands HME) for vocal and pulmonary rehabilitation after total laryngectomy. *Acta Otolaryngologica*. 2003 Jan;123(1):91-9.
- van As CJ, Hilgers FJM, Koopmans-van Beinum FJ, Ackerstaff AH. The influence of stoma occlusion on aspects of alaryngeal speech. *Acta Otolaryngologica*. 1998 Sep;118(5):732-8.

The above references represent a selection. For more information contact Atos Medical AB or your local representative.

Contact

Below you will find all our global addresses and contact information. Find your country and reach out to us. We are always here to help you.



Head office & Manufacturer:

Atos Medical AB
Tel. +46 (0)415 198 00
info@atosmedical.com
www.atosmedical.com

Australia:

Atos Medical Pty Ltd
Tel. +61 1800 286 728
info.au@atosmedical.com
www.atosmedical.com.au

Belgium:

Atos Medical BVBA / SPRL
Tel. +32 (0)2 218 55 50
info.be@atosmedical.com
www.atosmedical.be

Brasil:

Atos Medical ltd
Tel. +55 11 2305 2022
info.br@atosmedical.com
www.atosmedical.com.br

Canada:

Atos Medical Canada Inc
Tel. +1 833 514 2867
info.ca@atosmedical.com
www.atosmedical.ca

France:

Atos Medical S.A.S.
Tel. +33 (0)1 49 08 01 00
info.fr@atosmedical.com
www.atosmedical.fr

Germany:

Atos Medical GmbH
Tel. +49 (0)800 53 53 667
info.de@atosmedical.com
www.atosmedical.de

Hungary:

Atos Medical AB Magyarorszag
Tel. + 36 70 377 2414
www.atosmedical.com

Italy:

Atos Medical Srl
Tel. +39 049 775 522
info.it@atosmedical.com
www.atosmedical.it

Japan:

Atos Medical Japan Inc.
Tel. +81 (0)3 4589 2830
info.jp@atosmedical.com
www.atosmedical.jp

The Netherlands:

Atos Medical B.V.
Tel. +31 (0)79 593 5000
info.nl@atosmedical.com
www.atosmedical.nl

New Zealand:

Atos Medical Ltd
Tel. +64 (0)4 576 2100
info.nz@atosmedical.com
www.atosmedical.co.nz

Nordic countries:

Tel. +46 (0)415-19820
kundtjanst@atosmedical.com
www.atosmedical.com

Poland:

Atos Medical Poland Sp. z o.o.
Tel. +48 22 295 00 16
biuro@atosmedical.com
www.atosmedical.com

Portugal:

Atos Medical Spain S.L.
Sucursal em Portugal
Tel. +351 212 539 281
info.pt@atosmedical.com
www.atosmedical.pt

Spain:

Atos Medical Spain S.L.
Tel. +34 93 323 91 95
info.es@atosmedical.com
www.atosmedical.es

Switzerland:

Atos Medical, Switzerland
Tel. +41 (0)44 380 60 90
info.ch@atosmedical.com
www.atosmedical.ch

U.K.:

Atos Medical UK
Tel. +44 (0)115 784 1899
info.uk@atosmedical.com
www.atosmedical.co.uk

U.S.:

Atos Medical Inc.
Tel. +1 800 217 0025
info.us@atosmedical.com
www.atosmedical.us

**“When I got my first
Provox voice prosthesis
I could make myself
heard again.
It was an enormous
relief.”**

- Anna



At Atos Medical, we are committed to giving a voice to people who breathe through a stoma, with design solutions and technologies built on decades of experience and deep understanding of our users.

Head office:
Atos Medical AB, Hyllie Boulevard 17, SE-215 32 Malmö, Sweden
Tel: +46 (0)415 198 00, Fax: +46 (0)415 198 98
Email: info@atosmedical.com, Web: www.atosmedical.com

Manufacturer:
Atos Medical AB, Kraftgatan 8, P.O. Box 183, SE-242 22 Hörby, Sweden
Tel: +46 (0)415 198 00, Fax: +46 (0)415 198 98
Email: info@atosmedical.com, Web: www.atosmedical.com

© Atos Medical AB, Sweden.

Caution:
This catalog is intended for distribution outside the U.S. only.

A decorative graphic consisting of multiple overlapping, wavy red lines that create a sense of motion and depth, resembling a stylized sound wave or a series of ripples. It spans across the bottom half of the page.

Atos
Atos Medical **Your voice**